



[Procedure]

- (1) A dry, glass flask equipped with a Teflon-coated magnetic stirring bar and a condenser is evacuated and filled with an inert gas (argon or nitrogen) three times.
- (2) ZnTAC24[®] (23.8 mg, 0.0250 mmol calcd. as Zn₄(OCOCF₃)₆(O)) are placed into the flask under an inert gas.
- (3) The alcohol (2.0 mmol) and ethyl acetate (3.4 mL, 0.6 M) are transferred into the flask by syringe under an inert gas (Figure 1).
- (4) The reaction mixture is stirred under reflux condition for 18 h.
- (5) The reaction mixture is cooled to ambient temperature and concentrated. The residue is purified by silica gel column chromatography to give the desired acetate in 84% yield.




Figure 1.

[Notes]

Guaranteed grade solvent is distilled under an inert gas (argon or nitrogen atmosphere).

[Waste Disposal Information]

All toxic materials should be disposed of in accordance with "Prudent Practices in the Laboratory: Handling and Management of Chemical Hazards, Updated Version "; National Academy Press; Washington, DC, 2011. doi. [10.17226/12654](https://doi.org/10.17226/12654)

Visit Our Web Site !!!

2019/05/01



finechemicals@takasago.com

(C) 2019 TAKASAGO INTERNATIONAL CORPORATION - ALL RIGHTS RESERVED.